


# Ampli1™ KITS

## SINGLE CELL PRECISION



### Ampli1™ WORKFLOW: SINGLE CELL DNA ANALYSIS


Single-cell analysis allows the study of cell-to-cell variation within a cell population and is a prerequisite for in-depth understanding of tumor heterogeneity and to study diseases and drug development. Successful analysis of single cell genome starts with an appropriate Whole Genome Amplification (WGA) step and ends with an optimized sequencing process and bioinformatic analysis matching the characteristics of the adopted WGA. **Ampli1™ kits** provide a complete range of products covering the entire workflow and optimized for different kinds of genetic analysis on the amplified DNA.



### Ampli1™ WGA PLUS KIT: BEST-IN-CLASS FOR SINGLE CELL WHOLE GENOME AMPLIFICATION

The **Ampli1™ Whole Genome Amplification (WGA Plus) kit** has been developed and optimized for balanced and complete amplification of the genomic DNA content of a single cell.

Following cell lysis, DNA is **digested with a restriction enzyme** and adaptors are ligated onto the DNA fragments. Amplification is mediated by a **single highly specific PCR primer**.



## Ampli1™ WGA METHOD: LOW ALLELIC DROPOUT (ADO) RATE


Most WGA methods use random priming and uncontrolled amplification, causing frequent allelic dropout due to preferential (unbalanced) amplification of one of the two alleles.

**Ampli1™ WGA** provides a robust and reproducible amplification efficiency on every locus and on each one of the two alleles reducing preferential amplification and achieving the **lower ADO rate** compared to other methods.

%ADO Measured					
Papers WGA	2014 Blinder V et al WGS	2015 Huang L et al WGS	2016 Babayan A et al WES	2017 Borgstrom E et al WES	2016 Normand E et al STR
<b>Ampli1™</b>	<b>2%</b>		<b>9%</b>	<b>7-9%</b>	<b>6%</b>
<b>Competitor1</b>			24%	65-98%	49%
<b>Competitor2</b>		21-28%		18-47%	
<b>Competitor3</b>		76%			69%
<b>Competitor4</b>		33-38%	100%	93-95%	

*Blinder V et al, Hum Mutat. 2014 - Huang L et al, Annu Rev Genomics Hum Genet. 2015 - Babayan A et al, Oncotarget. 2016 - Borgström E et al, PLoS One. 2017 - Normand E et al, Prenat Diagn. 2016*

## Ampli1™ WGA METHOD: ROBUST AMPLIFICATION FROM BOTH FRESH AND FIXED CELLS



Ampli1™ WGA Method performs best **both with fresh and fixed cells** in terms of absolute DLRS value and of standard deviation compared to other WGA methods. Derivative Log Ratio Spread (DLRS) are indicative of the noise of the system. **Lower values allow higher resolution analysis of Copy Number Variation (CNV).**

*Normand E et al, Prenat Diagn. 2016, "Comparison of three whole genome amplification methods for detection of genomic aberrations in single cells"*

## NEW Ampli1™ WGA PLUS FEATURES:



**LOWEST ALLELIC DROP OUT**



**BEST IN CLASS FOR SINGLE CELL WGA**



**COMPATIBLE WITH LIQUID HANDLERS**



**NEW 4 HOURS PROTOCOL!**



**REPRODUCIBLE AMPLIFICATION EFFICIENCY**



**USER FRIENDLY PROCEDURE**

# Ampli1™ WGA QC KIT: QUALITY CONTROL OF WGA PRODUCT


The Ampli1™ WGA QC Kit provides a useful quality control check of the WGA procedure by utilizing a PCR-based assay to establish DNA integrity. Cells with good quality DNA typically produce 3 or 4 PCR bands, while cells with degraded DNA will show fewer QC bands. The number of bands defines a "Genome Integrity Index" (GII) for predicting success of various downstream analytical methods. (Polzer B et al., 2014 EMBO Molecular Medicine)

AMPLI1™ KITS ARE FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES  
The performance characteristics and safety and effectiveness have not been established and are not cleared or approved by the FDA.

# Ampli1™ LOWPASS: SPECIFIC, SENSITIVE AND QUANTITATIVE DETECTION OF COPY-NUMBER ABERRATIONS (CNAs) IN SINGLE CELL

The Library preparation kit is used to detect chromosomal aneuploidies and CNAs by lowpass whole genome sequencing. It is designed to work with the Ampli1™ WGA kit.

CK / CD45 / DAPI



Copy-Number Alteration patterns in cancer:

- Prognostic factors
- Markers of immune evasion
- Predictive factors of chemoresistance

Ferrarini et. al, PLOS One 2018 <https://doi.org/10.1371/journal.pone.0193689>

For Illumina and IonTorrent Platforms

# Ampli1™ ONCOSEEK PANEL: ACCURATE DETECTION OF SOMATIC VARIANTS AND FOCAL CNAs FROM SINGLE CELL

The Ampli1™ OncoSeek Panel enables the preparation of high-quality targeted NGS libraries for the simultaneous detection of single nucleotide variants (SNVs), insertions/deletions (indels) and copy-number amplifications (CNAs) in a panel of 60 clinically relevant oncology-related genes. It is designed to work with the Ampli1™ WGA product.

ABL1	CDK6	FBXW7	HRAS	MLH1	PTPN11	
AKT1	CDKN2A	FGFR1	IDH1	MPL	RG1	
ALK	CSF1R	FGFR2	IDH2	MSH6	RET	
APC	CTNNB1	FGFR3	JAK2	MYC	SMAD4	<ul style="list-style-type: none"> <li>• First and only panel designed for single cell analysis.</li> <li>• Coverage of 60 oncology-relevant genes including more than 2500 mutation hotspots and CNAs for a subset of 19 genes.</li> <li>• Very Low Allelic Dropout Rate (ADO=12.7% ± 4.2%)</li> </ul>
AR	DDR2	FGFR4	JAK3	MYCN	SMARCB1	
ATM	DNMT3A	FLT3	KDR	NOTCH1	SMO	
BRAF	EGFR	GNA11	KIT	NRAS	SRC	<p>For Illumina Platforms</p>
CCND1	ERBB2	GNAQ	KRAS	PDGFRA	TP53	
CDH1	ERBB4	GNAS	MAP2K1	PIK3CA	TSC1	
CDK4	EZH2	HNF1A	MET	PTEN	HVL	

# Ampli1™ WORKFLOW: CNV AND TARGET SEQUENCING FROM THE SAME SINGLE CELL

## Mutational and copy number status of APC gene in a patient with metastatic cancer of unknown primary

In CTC1 and CTC2 a APC mutation is detected by Ampli1™ OncoSeek in homozygous state (Fig1), which is consistent with the loss of one of the 2 alleles detected by Ampli1™ LowPass (Fig2).

The analysis shows a deletion of APC gene in CTC3 consistently detected by both Ampli1™ LowPass and OncoSeek (Laprovitera et al., 2021 Front. Cell Dev. Biol.):

- Ampli1™ LowPass detected a genome segment of chromosome 5 with Copy Number = 0
- Ampli1™ OncoSeek shows no coverage on all APC gene amplicons in agreement with the presence of a deletion.

gene	effect	hgvsp	CTC1	CTC2	CTC3
APC	frameshft_variant	p.T1556Nfs*3	90,22	99,23	0,00

Fig1: Ampli1™ OncoSeek somatic mutation detection in 3 CTCs of a patient with metastatic cancer of unknown primary. An APC mutation is detected only in CTC1 and CTC2.




Fig2: Ampli1™ LowPass genome-wide copy number profiles of CTC1, CTC2, CTC3 from the same patient represented in Fig.1

Ampli1™



MENARINI  
silicon biosystems

MENARINI SILICON BIOSYSTEMS, INC.  
3401 Masons Mill Rd Suite 100  
Huntingdon Valley, PA 19006 UNITED STATES  
Telephone: +1 877-837-4339

©2023 Menarini Silicon Biosystems S.p.A. All rights reserved  
Menarini Silicon Biosystems Spa | Via Giuseppe di Vittorio,  
21 b/3 40013 Castel Maggiore (BO) ITALY  
telephone +39 051 9944100 | fax +39 051 7459550

AMPLI1™ KITS ARE FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES  
The performance characteristics and safety and effectiveness have not been established and are not cleared or approved by the FDA.  
MKT-000096332